
Clement
Clarke
International

Flutter

Flutter

Instructions for use


Fig. 1



Fig. 2

Instructions for Use

Contents

	Page
Exploded Diagram Explanation	4
Intended Purpose	4
Intended User	4
Contraindications	4
How Flutter Works	4
Cautions	4
Preparation	5
Technique	6
Stage 1 - Loosening and motion	6
Stage 2 - Mucus removal with cough	7
Session Length and Frequency	7
Importance of Holding Breath	7
Cleaning Instructions	8
Life Expectancy of Flutter	8
Disposal Instructions	8
Symbols	87

Exploded Diagram Explanation

(Fig. 1)

1. a mouthpiece (stem)
2. a protective cover
3. a high-density stainless-steel ball
4. a circular cone

Intended Purpose

For use in the removal of excessive secretions from the lungs.

Intended User

For patients with cystic fibrosis, bronchitis and bronchiectasis, age 3+.

Contraindications

 Flutter should **not** be used by patients with the following conditions:

- Pneumothorax
- Right-sided Heart Failure
- Severe Tuberculosis
- Haemoptysis (coughing-up blood)

Consult your health care professional before use.

How Flutter works

When breathing out through the Flutter, the motion will cause the steel ball to bounce under the protective cover. The vibrations will loosen mucus sticking to your airways and begin to exit the lungs.

Cautions

- Not suitable for multi-patient use; due to risk of infection.
- Contains small parts. The Flutter is not a toy.
- If you have been taught a specific technique already, then continue to use it, as it will have been chosen as the best method for you.
- Regularly check that your technique is correct with your healthcare professional or by re-reading these instructions.
- Any serious incident in relation to this device must be reported to the manufacturer and the Competent Authority of the country where the event occurred.

Preparation

Ensure device is correctly assembled

1. You may want a glass of water and box of tissues within reach.
2. Begin by relaxing.
 - Sit with your back straight and head tilted slightly upwards so your throat and windpipe are wide open. This allows the air you breathe out to flow smoothly from your lungs out through the Flutter. Or you may prefer an alternative position, with elbows placed on a table at a height that won't allow you to 'slouch'. Experiment to find out which position is best (see Fig. 2).
 - If you must perform Flutter therapy while in bed, sit as upright as possible in a position that will not restrict your smooth breathing effort.
3. Keep your cheeks as stiff as possible, since the vibrations produced by the Flutter are wasted when you allow your cheeks to move a lot. Look into a mirror or hold your cheeks flat with your spare hand to help learn the technique. Eventually you will learn what to look out for and your cheeks will not vibrate as much.
4. The angle at which you hold the Flutter is critical. Start by holding it level, with the stem horizontal to the floor. This places the cone at a slight tilt and insures that the ball not only flutters but also rolls during exhalation. Experiment by tilting the Flutter upward and downward slightly (see Fig. 3). Use the angle that gives you the strongest fluttering sensation. The fluttering may not occur immediately.

Technique

Stage 1 - Loosening and motion

1. Slowly, breathe in beyond a normal breath, do not fill your lungs completely.
2. Hold your breath for 2 to 3 seconds.
3. Place the Flutter in the mouth with lips closed around the stem, making sure no air can escape around the mouthpiece.

Do not cover the exit holes while holding the Flutter.

4. Position the Flutter at the desired angle.
5. Begin to breathe out through the Flutter at a reasonably fast but not too forceful speed.

Exhaling too fast and too forcefully at this point should be avoided.

Concentrate on feeling the vibration in your airways, since the best results are obtained when there is maximum fluttering in your chest.

6. Continue to exhale through the Flutter to a level slightly further than you would exhale during normal breathing.

For the moment, any urge to cough should be suppressed.

7. Repeat this breathing technique for 5 to 10 breaths to loosen and move as much mucus as possible from your airways. The mucus will move further up the airways with each Flutter exhalation. It is important to inhale slowly (away from the Flutter or through your nose) and hold your breath, for 2 to 3 seconds, before each exhalation through the Flutter.

For the moment, any urge to cough should still be suppressed.

Inhaling slowly helps to avoid mucus being sucked back into the lungs.

'First-time' users or those with particularly thick mucus may find that repeating stage 1 many times is necessary before starting stage 2. Several attempts may be necessary to help determine the appropriate number of breaths you need to perform at each stage.

Stage 2 - Mucus removal with cough

1. Perform 1 or 2 additional breaths through the Flutter.
2. Breathe in as slowly and fully as comfortable, filling your lungs with air.
3. Hold your breath for 2 to 3 seconds.
4. Breathe out through the Flutter as forcefully and fully as you comfortably can.

This forceful exhalation moves mucus up to a level in your lungs that triggers a cough. If mucus is not easily coughed out following 1 or 2 attempts, try a 'huff' manoeuvre (like trying to remove a bread crumb from your throat).

Usually, you will need to do 5 to 10 mucus loosening and motion breaths, suppressing a cough, followed by 1 or 2 mucus removal breaths, to clear your airways successfully.

If you find your requirements for this therapy (loosening, moving and removing mucus) changes, you may need to modify your routine to reach the best therapeutic effect.

Session Length and Frequency

Two sessions a day (morning and evening) are recommended. Sessions will usually range from 5 to 15 minutes, depending on the condition of your lungs.

To avoid becoming too tired, you can take three sessions instead of two per day. The inability to cough out additional mucus during therapy often indicates the session has been successful and can be ended.

Importance of Holding Breath

People with lung conditions cannot fill their lungs evenly. Holding your breath for 2 to 3 seconds allows the inhaled air to be distributed more evenly throughout the lungs.

Cleaning Instructions

Wash the Flutter at least every 2 days.

Unscrew protective cover and separate components (see Fig. 1).

Your Flutter device can be cleaned by immersing the device in warm (but not hot) mild detergent solution for 2-3 minutes (maximum 5 minutes). Agitate the device to ensure thorough cleaning. Rinse in warm water and shake to remove any excess water. It is important to rinse thoroughly to prevent corrosion of the stainless-steel ball. Dry with a clean towel. When dry, reassemble and store in a clean, dry location.

Ensure the device is thoroughly dry before storage.

Do not clean in a dishwasher. This will cause damage to the components of the device. (e.g. rust).

Do not use chlorine, bleach or other chlorine-containing products (e.g. dishwasher machine detergent). This will cause damage to the components of the device (e.g. rust).

Life Expectancy of Flutter

The Flutter should be replaced every 12 months or sooner if damaged.

Disposal Instructions

Follow cleaning instructions and dispose in household waste.

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Made by Clement Clarke International Ltd. in Mountain Ash, United Kingdom, under licence from ACTAVIS Group PTC EHF. FLUTTER is a registered trademark of ACTAVIS Group PTC EHF.
Clement Clarke International Ltd.
Flexicare Deutschland GmbH
Flexicare Medical Switzerland GmbH
c/o BK-Services AG Suurstoffi 22 • 6343 Rotkreuz Switzerland
Clement Clarke International Ltd is part of the Flexicare Group.