### Instructions for Use—English

#### Model 8000AP Pediatric Finger Clip Pulse Oximeter Sensor

**Indications for Use**  
The Nonin Model 8000AP Pediatric Finger Clip Sensor is designed for spot-to-spot monitoring of oxygen saturation monitoring of pediatric patients where little sensor motion is expected.

**CAUTION:** Federal law (USA) restricts this device to sale by or on the order of a licensed practitioner.

#### Warnings:

- Do not use the device in an MR environment.
- Do not use on Nonin pulse oximeters. These pulse oximeters are manufactured to meet the accuracy specifications for Nonin sensors. Using other pulse oximeters may cause improper sensor performance.
- Inspect the sensor application site at least every 4 hours to ensure correct sensor alignment and skin integrity. Patient sensitivity to sensors may vary due to medical issues or skin condition.
- Avoid excessive pressure to the sensor application site as this may cause damage to the skin beneath the sensor.

- ⚠ Cautions:**
- Do not prevent improper performance and/or patient injury, verify sensor and pulse oximeter compatibility before use.
  - Do not use a damaged sensor. If the sensor is damaged in any way, do not use, or discontinue use immediately.
  - Follow local governing ordinances and agency instructions regarding disposal or return of the sensor and any components.
  - A functional tester cannot be used to assess the accuracy of a pulse oximeter of another or probe.
  - Refer to the pulse oximeter operator's manual for additional warnings and cautions.
  - Factors that may degrade pulse oximeter performance include the following:
    - excessive ambient light
    - anemia, polycythemia
    - radioactive interference
    - relaxed catheters, blood pressure cuffs, infusion lines
    - moisture in the sensor
    - improperly applied sensor
    - carboxyhemoglobin
    - methemoglobin
    - artificial nails
    - incorrect sensor type

Symbols:	Symbol	Definition	Symbol	Definition
		Follow instructions For Use.		Indicaciones separadas de limpieza para el dispositivo electrónico y el equipo médico.
		CAUTION		Storage/preparation temperature range
CE 0123		CE Marking indicating conformance to IEC Directive No. 60427/IEC concerning medical devices		Manufacturer
		Protected against vertically falling water drops when enclosure is tilted up to 15 degrees and ingress of solid foreign objects larger than or equal to 2.5 mm in diameter per IEC 60529.		Country of origin
IP32		IP32		Date of manufacture
		Lot		Catalogue number
		RoHS compliant (China)		Quantity
		Authorized representative in the European Community		RoHS compliant (China)
		Unique Device Identifier		Medical prescription required
		Medical Device		Keep Dry
		Distributor		Importer
		Use By		Handle With Care

### Choosing the Appropriate Sensor

The 8000AP is designed for use on the fingers of patients weighing between 8 and 30 kilograms, where the finger tissue thickness is between 5 and 13 millimeters (dimension "H" in Figure A). An excessive sensor motion is expected with the Model 8000A Infant Finger Sensor.

### Attaching the Pediatric Finger Clip Sensor

1. Insert a finger (preferably the index, middle, or ring finger) into the Pediatric Finger Clip Sensor (Figure B) until the end of the finger reaches the finger stop. Keep the finger nail facing the sensor (as shown in Figure B). Ensure that the finger nail does not interfere with proper finger position.

2. Fasten the sensor to the finger by pulling the sensor strap across the monitoring, secure the sensor cable independently from the sensor with medical tape, preferably around the base of the finger. Make sure that the tape securing the cable does not restrict the blood flow.

The thumb is not recommended for use with the Pediatric Finger Clip Sensor.

**Note:** Sensor placement is critical for good performance. If the sensor is not positioned properly, light may bypass the tissue and result in SpO<sub>2</sub> inaccuracies.

**Cleaning the Reusable Sensor**  
**⚠ Cautions:**

- Do not clean the sensor before applying it to a new patient.
- Unplug the sensor from the pulse oximeter before cleaning.
- Do not sterilize, autoclave, or immerse the sensors in liquid of any kind. Do not spray any liquids into the sensor. Do not sterilize with EO.
- Do not use caustic or abrasive cleaning agents on the sensors. Do not use cleaning agents containing ammonium chloride or isopropyl alcohol.

1. To do in cleaning, remove the spring from the slot on the bottom of the sensor (Figure C), and then rotate the spring forward as shown in Figure D.

**Note:** The spring must be disengaged from the bottom slot before opening the sensor 90° or the spring will be damaged.

2. Clean the sensor, wipe all patient contact surfaces with a soft cloth dampened with a mild detergent or a 10% bleach/90% water solution (household bleach [containing less than 10% sodium hypochlorite]). Ensure that all tape residues is removed.

3. Reverse the steps to thoroughly before reusing.

4. Reverse the steps shown in Figures C and D to reattach the sensor spring.

**Note:** To minimize cable deterioration when cleaning the cable, gently wipe away from the plug end towards the sensor end.

Specifications SpO <sub>2</sub> Accuracy 1, 2	Oxygen Saturation (A <sub>max</sub> ) (Figure 1)	Motion Oxygen Saturation (A <sub>min</sub> ) (Figure 2)
70 - 100%	±2	±3
70 - 80%	±2	±3
80 - 90%	±2	±3
90 - 100%	±2	±2

**SpO<sub>2</sub> Low Perfusion Accuracy:** 70% to 100% ±2 digits (A<sub>max</sub>)  
**Pulse Rate Accuracy:** 40 to 240 BPM ±3 digits (A<sub>max</sub>)  
**Pulse Rate Motion Accuracy:** 40 to 250 BPM ±3 digits (A<sub>max</sub>)  
**Pulse Rate Low Perfusion Accuracy:** 40 to 240 BPM ±3 digits (A<sub>max</sub>)

**Temperature:** 3, 4  
Operating: -20 °C to 50 °C (-4 °F to 122 °F)  
Storage/Transportation: -40 °C to 70 °C (-40 °F to 158 °F)  
Humidity: 3, 4  
Operating: 10% to 95% non-condensing  
Storage/Transportation: 10% to 95% non-condensing

**RoHS:** 860 nanometers @ 3 mW nominal  
Infrared: 910 nanometers @ 3 mW nominal  
*These information is especially useful for clinicians.*

**Compliance**  
This product complies with ISO 10993.  
Not made with natural rubber latex.

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Users and/or patients should report adverse events involving their Nonin device to Nonin Medical, Inc. and the competent authority of the EU Member State in which the user and/or patient is established, if applicable.

**Warranty**  
1 year from the date of delivery.  
The device's expected service life is 1 year.  
Nonin reserves the right to make changes and improvements to these instructions and the product it describes at any time, without notice or obligation.  
www.nonin.com/warranty

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