

# Stop voor twee



Baby CO



CO (ppm)	Moeder	%COHB	%FCOHB	Foetus
20 & above		3.20+		6.26+
20		3.20		6.26
19		3.04		5.95
18		2.88		5.63
17		2.72		5.32
16	<b>Roker</b>	2.56		5.00
15		2.40		4.69
14		2.24		4.38
13		2.08		4.07
12		1.92		3.76
11		1.76		3.44
10		1.60		3.13
9	<b>lichte roker</b>	1.44		2.82
8	<b>roker</b>	1.28		2.50
7		1.12		2.19
6		0.96		1.88
5	<b>Niet roker</b>	0.80		1.56
4	<b>roker</b>	0.64		1.20
3		0.48		0.94
2		0.32		0.63
1		0.16		0.31

**Hoge risico zone**

**Gezond**