



Medical

CO (ppm) **%COHB**

20 en hoger 3.20+

20 3.20

19 3.04

18 2.88

17 2.72

16 2.56

15 2.40

14 2.24

13 2.08

12 1.92

11 1.76

10 1.60

9 1.44

8 1.28

7 1.12

6 0.96

5 0.80

4 0.64

3 0.48

2 0.32

1 0.16

Alarm

Zware roker

Lichte roker

Niet roker

De Smoke Check en Micro CO koolmonoxidemonitoren en accessoires maken deel uit van een uitgebreide reeks respiratoire meetapparatuur.

Ze zijn te bestellen aan de hand van onderstaande omschrijvingen:

Product

Micro CO:	MC02
Smoke Check:	SC01
Micro CO kalibratiegas:	MCG020
Smoke Check kalibratiegas:	MCG020
Kalibratieventiel met stromingsindicatie:	MGA222
COBRA softwarepakket inclusief aansluitkabel voor Micro CO:	CBI1000
Nieuwe sensor voor Micro CO:	CEL5900
Interne lithiumbatterij voor Micro CO:	BAT5900
Nieuwe sensor voor Smoke Check:	CEL4100
Mondstukadapter met eenwegklep	PSA2000
Kartonnen mondstukken, 22mm voor eenmalig gebruik (100 per doos)	
Kartonnen mondstukken, 22mm voor eenmalig gebruik (500 per doos)	

Micro Medical Ltd pursues a policy of continuing improvement in design, production and performance of its products, the right is therefore reserved to vary details at any time without notice.



PT Medical BV
Lorentzpark 12, 9351 VJ Leek
Tel 0594-587 280
Fax 0594-587 288
Email mail@pt-medical.nl
www.pt-medical.nl



Micro Medical Limited
PO BOX 6, Rochester, Kent, ME1 2AZ, UK
Tel +44 1634 893500
Fax 01634 893600
Email micromedical@viasyshc.com
www.micromedical.co.uk

