

# Stop voor twee

CO Check Baby



| CO (ppm)   | Moeder              | %COHB | Foetus |
|------------|---------------------|-------|--------|
| 20 & above |                     | 3.20+ | 6.26+  |
| 20         |                     | 3.20  | 6.26   |
| 19         |                     | 3.04  | 5.95   |
| 18         |                     | 2.88  | 5.63   |
| 17         |                     | 2.72  | 5.32   |
| 16         | <b>Roker</b>        | 2.56  | 5.00   |
| 15         |                     | 2.40  | 4.69   |
| 14         |                     | 2.24  | 4.38   |
| 13         |                     | 2.08  | 4.07   |
| 12         |                     | 1.92  | 3.76   |
| 11         |                     | 1.76  | 3.44   |
| 10         |                     | 1.60  | 3.13   |
| 9          | <b>lichte roker</b> | 1.44  | 2.82   |
| 8          | <b>roker</b>        | 1.28  | 2.50   |
| 7          |                     | 1.12  | 2.19   |
| 6          |                     | 0.96  | 1.88   |
| 5          | <b>Niet roker</b>   | 0.80  | 1.56   |
| 4          | <b>roker</b>        | 0.64  | 1.20   |
| 3          |                     | 0.48  | 0.94   |
| 2          |                     | 0.32  | 0.63   |
| 1          |                     | 0.16  | 0.31   |

**Hoge risico zone**

**Gezond**